

## Test Report

Prepared for:

Lance Liberti, MBA  
Chief Executive Officer  
Integrative Practice Solutions

Dr. Robert McGrath DO  
Medical Director  
Integrative Practice Solutions  
Juventix Regenerative Medical

01 March 2023

PN:2301

**I. Objective of Study:** The primary objective of this study was to evaluate platelet rich fibrin (PRF) product following an activation process. A secondary objective was to measure pH and platelet concentration of platelet rich plasma (PRP) products.

### III. Study Parameters

#### *Platelet Concentration*

Complete blood counts (CBCs) were performed using a hematology analyzer for the baseline and concentrate samples.

#### *pH*

pH was measured for the PRP sample using a blood gas analyzer.

#### *Platelet Derived Growth Factor Analysis*

Growth Factor quantitation was performed by ELISA for untreated and activated PRP samples. A positive control sample prepared using CaCl<sub>2</sub> and bovine thrombin was also tested.

**Table 3. Summary of Spin Protocols for 15mL PRP Tube**

Spin Protocol	Centrifugation Speed, Time	Collection Method
Protocol	3500 rpm, 8 minutes	immediate removal of all but ~2 mL of plasma post-centrifugation from top of tube; 20-minute rest period with remaining 2mL; gentle mixing to resuspend platelets then collection of 2mL of PRP

**Table 4. Hematology Data – 15 mL Tube**

		Protocol	
Cell Count Parameters		Whole Blood	PRP
WBC	(x 10 <sup>6</sup> /mL)	3.4	0.75
RBC	(x 10 <sup>9</sup> /mL)	4.0	0.03
PLT	(x 10 <sup>6</sup> /mL)	161	548
LY	(x 10 <sup>6</sup> /mL)	1.42	0.66
MO	(x 10 <sup>6</sup> /mL)	0.26	0.04
NE	(x 10 <sup>6</sup> /mL)	1.62	0.03
EO	(x 10 <sup>6</sup> /mL)	0.10	0.01
BA	(x 10 <sup>6</sup> /mL)	0.02	0.00
<b>PLT Recovery (%)</b>			<b>84%</b>

**pH**

**7.5**

*WBC – White Blood Cell; RBC – Red Blood Cell; PLT – Platelet; LY – Lymphocyte; MO – Monocyte; NE – Neutrophil; EO – Eosinophil; BA – Basophil*

**Table 5. Platelet Derived Growth Factor Concentration (pg/mL)**

<b>Test Samples</b>	<b>PDGF (pg/mL)</b>
Untreated PRP	1536
Treated PRP	10844

## **V. Summary of Results**

This study assessed parameters associated with PRP and PRF products prepared using Juventix PRP devices and the photoactivation/bio-incubation method. PRP products were prepared from 15mL and 30mL PRP tubes using several centrifugation protocols and platelet concentrate collection methods. Substantial platelet recoveries (84% - Table 4) were calculated for the 15mL PRP device, when whole blood was processed using centrifugation protocol, respectively. Both protocols included a 20-minute 'rest' period, which has been previously shown to improve platelet yields (*PN2006, July 2020 study*). The measured pH of platelet concentrates was ~7.5.

PDGF growth factor concentrations were determined for untreated PRP and photoactivated/bio-incubated PRP/PRF test samples. A CaCl<sub>2</sub>-Thrombin generated releasate prepared from the same PRP product was included in the analysis as a positive control. Photoactivation and bio-incubation led to the transformation of the liquid PRP product to a gel-like PRF product, which indicates platelet activation occurred and the clotting cascade was initiated.